

true training
true anatomy
true to life

AIRSIM ADVANCE 

tru
corp



TruCorp Ltd
The Mount Business Park
2 Woodstock Link
Belfast BT6 8DD

T. +44 028 9073 7281
F. +44 028 9073 7282
E. info@trucorp.co.uk

www.trucorp.co.uk



AIRSIM ADVANCE

The AirSim Advance features the uniquely constructed AirSim airway, designed to provide true, anatomically correct and visually accurate internal features.

The addition of a “real-feel” skin covering provides a more realistic bag mask ventilation training experience. AirSim Advance features an improved neck design allowing for accurate articulation.

PRODUCT FEATURES

AirSim Advance also features an anatomically accurate nasal passage, inflatable tongue to increase difficulty of airway and an option to include “breakout teeth” to simulate the effects of bad practice in laryngoscopy.

AirSim Advance is ideal for training in:

- Effective bag mask ventilation.
- Use of full range of supraglottic devices.
- Video and standard laryngoscopy.
- Tracheal & naso-tracheal intubation.
- Awake fiberoptic examination.

Trucorp recommend the following equipment sizes for optimal performance:

- 7.0-7.5 mm ID for nasal intubation
- 8.0-9.0 mm ID for oral intubation
- Size 3 and 4 for LMA laryngeal masks (similar respective sizes for other supraglottic devices)
- Size 35F - 37F for endobronchial tubes.

Product Name	Code	Description
AirSim Advance	AA91100	AirSim Advance Model
AirSim Advance	AA91101	Advance, Breakout teeth
AirSim Advance	AA95100	Advance with Bronchial Tree
AirSim Advance	AA95101	Advance, Bronchial Tree, Breakout teeth

