


true training
true anatomy
true to life

tru
corp

AIRSIM ADVANCE BRONCHI 



TruCorp Ltd
The Mount Business Park
2 Woodstock Link
Belfast BT6 8DD

T. +44 028 9073 7281
F. +44 028 9073 7282
E. info@trucorp.co.uk

www.trucorp.co.uk



AIRSIM ADVANCE BRONCHI

The AirSim Advance Bronchi is the next generation to the AirSim family, providing anatomically correct detail down to the fourth generation bronchi together with the addition of our 'real feel' skin covering. This model combines together the skin covering of the AirSim Advance with all the combined features of the AirSim Bronchi to provide a true, anatomically correct and visually accurate Bronchi trainer. The AirSim Advance Bronchi features an improved neck design allowing for accurate articulation and provides exceptional detail in both internal and external features

PRODUCT FEATURES

The AirSim Advance Bronchi can also be used in a similar manner to the AirSim Bronchi and can cover an extensive range of training requirements including Naso and Endo tracheal intubation.

The AirSim Advance Bronchi unit will provide a solution to your training needs in:

- Diagnostic bronchoscopy.
- Lung isolation techniques using left and right endobronchial tubes and bronchial blockers.
- Lung suctioning techniques.

For clinicians and industrial clients, the AirSim Advance Bronchi is a leading choice for demonstration of items such as:

- Fibreoptic bronchoscopic equipment.
- Single and double lumen tubes and bronchial blockers.
- Specialised bronchial stents.

Trucorp recommend the following equipment sizes for optimal performance:

7.0-7.5 mm ID for nasal intubation
8.0-9.0 mm ID for oral intubation

Size 3 and 4 for LMA laryngeal masks
(similar respective sizes for other supraglottic devices)
Size 35F - 37F for endobronchial tubes.

Product Name	Code	Description
AirSim Advance Bronchi	AA95100	AirSim Advance Bronchi model
AirSim Advance Bronchi	AA95101	AirSim Advance Bronchi with Breakout teeth

