

true training
true anatomy
true to life

tru
corp

AIRSIM 

ADVANCE
CRICO



TruCorp Ltd
The Mount Business Park
2 Woodstock Link
Belfast BT6 8DD

T. +44 028 9073 7281
F. +44 028 9073 7282
E. info@trucorp.co.uk
www.trucorp.co.uk



AIRSIM ADVANCE CRICO

The AirSim Advance Crico combines all the standard features of the Advance range with the addition of a simulated cricoid cartilage to facilitate training in both needle and surgical cricothyroidotomy procedure.

The AirSim Advance Crico is the solution for delivery of effective training in a comprehensive range of airway management procedures.

PRODUCT FEATURES

The AirSim Advance Crico features the uniquely constructed AirSim airway designed to provide true, anatomically correct and visually accurate internal features, coupled with a “real-feel” skin covering.

The AirSim Advance Crico is ideal for training in:

Needle/Surgical Cricothyroidotomy.

Awake fiberoptic examination.

Tracheal & naso-tracheal intubation.

Video and standard laryngoscopy - including Sellick Manoeuvre.

Use of full range of supraglottic devices.

Effective Bag Mask Ventilation.

Trucorp supply the following consumables to ensure optimal performance of the Advance Crico:

RS100X - Replacement neck skin attachments

LA1000 - Replacement larynx

ST1000 - Sealing tape for trachea

Product Name	Code	Description
Advance Crico	CR91100	AirSim Advance Crico model
Advance Crico	CR91101	AirSim Advance Crico model with Breakout teeth
Advance Crico	CR95100	AirSim Advance Crico model with Bronchial tree
Advance Crico	CR95101	AirSim Advance Crico model with Bronchial tree & Breakout teeth

