



# Parker Flex-It™ Directional Stylet





www.parkermedical.com



To Life!™

The Parker Flex-It™ Directional Stylet is an articulating stylet whose curvature may be selected and fixed **before** intubation, or controlled by easy thumb action **during** intubation. Its convenient, one-hand operation is designed for quick, accurate, gentle intubation. It is suitable for **oral** and **nasal** use.

#### Order #

-  FLEX-7580 For use in 7.5 and 8.0 tracheal tubes.
-  FLEX-6570 For use in 6.5 and 7.0 tracheal tubes.
-  FLEX-5560 For use in 5.5 and 6.0 tracheal tubes.
-  FLEX-5000 For use in 5.0 tracheal tubes.